

Fast + Simple
Focused on Veterinary Diagnostics

FASTest® CPV Ab ad us. vet. **FASTest® CDV Ab** ad us. vet.

Specific antibody and vaccination diagnostics

Fast test for the **qualitative** detection of antibodies against **Canine Parvovirus** or **Distempervirus** antibodies in whole blood, plasma or serum of the dog

Immune status of the breeding bitch
during pregnancy


**Determination of individual
vaccination time point**

- primary vaccination (basic immunisation)
- booster vaccination

Control of vaccination success

Control of vaccination titre
before travelling, exhibitions,
shelter dwelling etc.

Prognostic diagnostics
immune status in acute
CPV or CDV disease cases



■ Simple test procedure with whole blood, plasma or serum

■ Fast test interpretation after 10 minutes

■ Reliable clinical diagnostics

	Sensitivity	Specificity
CPV Ab	99.9%	94.0%
CDV Ab	99.9%	99.8%

■ Storage at room temperature (15-25°C)

■ Long shelf life

■ Compact test box with 2 or 10 tests



FASTest® CPV Ab_{ad us. vet.}

FASTest® CDV Ab_{ad us. vet.}

Antibodies are basic modules of the humoral immune response. They are passed by passively via the colostrum as so-called maternal antibodies (mAb) onto the yet immunoincompetent newborns or induced actively by natural field infection or vaccination. The antibody titre is varying individually in each animal, depending on multiple factors. The titre can persist over an extended period of time, partially lifelong, in efficient protection concentration (= reliable immunity by protective antibodies) or can fall below the efficient protection concentration (non-reliable immunity) in the course of time.

Each cut off (reliable immunity or not) of **FASTest® CPV Ab (1:80)** and **FASTest® CDV Ab (1:16)** is considered according to the Golden Standard Tests (haemagglutination inhibition test and virus neutralization test, respectively).

Depending on the level of individual antibody titre, the veterinarian is able to decide fast and reliable the necessity of vaccination or non-vaccination due to following questions:

Individual vaccination point

- of the breeding bitch

In problematic breedings, the determination of antibody status of the female makes sense during pregnancy to decide whether a booster vaccination before birth is necessary or to find the optimal primary vaccination time of the puppies.

- of the puppies: primary vaccination

There is a critical stage (so-called immunity gap) in puppies, especially in the first 12 weeks. During this stage the concentration of mAb could be high enough to inactivate the vaccinating virus but also too low to protect from field infection. Therefore it is important to find the individual primary vaccination point for each puppy to guarantee an appropriate protection.

For the determination of antibody status of the whole litter, it is possible to determine the antibody status of only one puppy, representative for the other puppies (so-called "fraternal antibody titre"). Here, the balanced colostrum assumption or development of all puppies is absolutely necessary.

- booster vaccination

By determination of the actual antibody status, an individual decision of the necessity of booster vaccination of the puppy or the adult animal can be made.

Being fast, safe and reliable, for pet owner and breeder these important questions can be answered practically by **FASTest® CPV Ab** or **FASTest® CDV Ab**. This enables the veterinarian an appropriate and customized vaccination diagnostics and strategy, adapted to dog and pet owner.

Test procedure

5 µl whole blood plasma serum

10 min

Test interpretation

ANTIBODY TITRE equal or higher than protective titre CPV Ab ≥ 1:80 CDV Ab ≥ 1:16	ANTIBODY TITRE lower than protective titre CPV Ab < 1:80 CDV Ab < 1:16
High titre colour intensity of T line > C line	Low titre colour intensity of T line < C line
Median titre colour intensity of T line ≈ C line	Titre below detection limit no T line
Good to very good CPV/CDV immune status, NO VACCINATION REQUIRED	Bad or no CPV/CDV immune status, VACCINATION RECOMMENDED

After positive antigen detection via **FASTest® PARVO Strip** and/or **FASTest® DISTEMPER Strip**, by application of **FASTest® CPV Ab** and/or **FASTest® CDV Ab** the antibody status of the diseased animal can be determined. This gives knowledge to the veterinarian about the immune status of the animal and therefore a better prognostic valuation of the further course of disease. And, what's more, the combined diagnostics alleviates the decision of therapy.

Distribution:

EN 03-2016

