

[Click Here For More Information About](#)

FASTest® IgG bovine

ad us. vet.



In vitro diagnosticum



Test-kit for the qualitative detection of IgG antibodies in whole blood, plasma or serum of cattle

INSTRUCTIONS FOR USE



Supplied Exclusively To The UK
Veterinary Market By
Vetlab Supplies Ltd
Visit Our Website
www.vetlabssupplies.co.uk
Telephone: 01798 874567
email us: info@vetlabssupplies.co.uk

Manufacturer:



1. INFORMATION ON THE TEST-KIT

TEST-KIT COMPONENTS

1 test-kit **FASTest® IgG** bovine contains:

- 2 or 10 test cassettes coated with monoclonal antibodies
- 1 dropper bottle **A** with 1.0 ml or 3.0 ml buffer diluent
- 2 or 10 disposable plastic pipettes
- 1 instructions for use

STABILITY AND STORAGE



Store at
15–25 °C



Expiry date
– see label

APPLICATION AND ABBREVIATIONS



For veterinary use only



Lot number



In vitro diagnosticum



Do not use test-kit components from different kits, lot numbers or beyond stated expiry date.



Follow instructions for use precisely

T – TEST line, **C** – CONTROL line, **LF** – Lateral flow

LIABILITY

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from the use of this product.

ACCURACY

Sensitivity 95.5%

Specificity 96.2%

(Comparison Method: Single Radialimmunodiffusion)

2. INTRODUCTION

An effective immunity is based on optimal levels of immunoglobulins (IgG, especially IgG1). This guarantees an optimal breeding success, especially in the peripartur stage and in new born calves, and is of economic importance to the livestock industry.

Due to the special placenta conditions, newborn calves often show no to hardly noteworthy IgG. Therefore the most important base for an immunoprophylaxis is an adequate supplementation with IgG.

IgG transfer mainly is based on the ingestion and resorption of colostrum (adequate passive transfer) directly after birth up to maximal 24h later. Failure of passive transfer (FPT) may be caused by inadequate suckling (reduced vitality, neonatal respiratory depression a.s.o.), way of colostrum application, very low levels of IgG in the colostrum, inadequate absorption of IgG or environmental stress. The consequences of an immune deficiency syndrome (IDS) are weak calves with a higher vulnerability for infectious newborn diseases like enzootic pneumonia of calves, diarrhoea (neonatal calf diarrhoea; calf scours) and other septicæmic diseases.

The incidence for cows (dystocia, metritis, mastitis a.s.o.) in the peripartur period is as high as in no other farm animal. Latest investigations show following IgG1 concentrations in healthy cows with normal parturition and physiologic postpartur period:

- Ø 38 mg/ml at 7th month of gestation
- Ø 15 mg/ml physiologically low concentration peripartur
- Ø 30 mg/ml increasing from 4th month post partum

If the intrapartur IgG concentration is already significantly decreased (< 15 mg/ml), an increased incidence of genital associated diseases will be seen in the peripartur period.

The optimal test time slot for using **FASTest® IgG** bovine is between 24 to 48h (max up to 7 days) post natum in calves and from 3rd to 7th day post partum in cows. **FASTest® IgG** bovine enables the veterinarian on-farm and without any technical costs to confirm (**FASTest® IgG** bovine: IgG ≤ 10 mg/ml) or to exclude (**FASTest® IgG** bovine: IgG > 10 mg/ml) a suspicion for FPT/IDS in the cow as well as in the calf.

3. INFORMATION ON THE SPECIMEN MATERIAL

Exactly 20 µl (of attached plastic pipette) 15–25 °C warm whole blood (WB, with anticoagulant), plasma (P) or serum (S) are needed.

Mix the sample material well before use!

Non-cooled (15–25 °C), WB, P and S should be tested immediately! At 2–8 °C, WB, P and S can be stored up to 4 days. Serum and/or plasma samples can be permanently stored at minimum –20 °C.

Keep in mind that the sample material, as well as all used test-kit components, should have reached **room temperature** at the time of application.

Endogeneous and exogeneous interfering substances of the sample (e.g. albumin, fibrinogen, lipids, CRP, heterophilic antibodies, especially type IgA, as well as viscosity, pH-value and excess EDTA) **as well as native blood can cause interferences** (matrix effects) **that can influence the target measurement. These can lead to an impaired LF and/or unspecific reactions on T and C.**

Whole blood samples have a lower sensitivity in the FASTest® IgG bovine. This can lead to false negative results!

4. SPECIMEN PREPARATION

- No specimen preparation necessary.
- **Partially filled and/or insufficient mixed EDTA, citrate or heparin tubes could create invisible microclots resulting in lateral flow delay and/or unspecific reactions** (e.g. greyish shadow like lines).

5. TEST PROCEDURE

1. Remove the test cassette from its foil pouch shortly before use. Place it on a flat surface.
2. Draw sample up to the mark (± 20 µl sample volume) using the disposable plastic pipette. **The meniscus must be above the black line** (fig.1).
3. Place the whole sample volume (20 µl) into the sample window **S** of the test cassette (hold pipette vertically, fig.2). Let the sample absorb completely – this can take up to a minute (whole blood)!
4. Hold the dropper bottle **A** vertically and express **3 drops of buffer solution** (ca. 120–150 µl) without bubbles into the sample window of the test cassette (fig.3).




fig.1




fig.2



fig.3

6. READING OF THE TEST RESULT



The test result must be read **exactly 10 minutes** after addition of the buffer into the sample window **S**. Beyond this time, test results should not be interpreted!

Negative: < 10 mg/ml ← cut-off → **Positive:** ≥ 10 mg/ml

NEGATIVE TEST RESULT: IgG < 10 mg/ml (fig.4)

A pink-purple **TEST line of any intensity** (varying from very weak to strongly intensive) and a pink-purple **CONTROL line** appear.

→ Suspicion on FPT/lack of IgG (IDS)

→ Colostrum administration is recommended

POSITIVE TEST RESULT: IgG ≥ 10 mg/ml (fig.5)

No pink-purple **TEST line** visible. Only the pink-purple **CONTROL line** appears. This line indicates, irrespective of its intensity, that the test has been performed properly.

→ adequate passive transfer

INVALID TEST RESULT

No **CONTROL line** visible. The test should be repeated using a new test cassette.

fig.4
NEGATIVE TEST RESULT



fig.5
POSITIVE TEST RESULT



7. PRECAUTIONS FOR USERS

- The guidelines for working in medical laboratories must be observed. It is recommended to wear disposable gloves and other personal protective equipment (protective clothing, possibly a face mask). Wash and disinfect hands after completing the test.
- Label sample material and associated test cassette to ensure a precise assignment.
- Use a new pipette and a new test cassette for each sample.
- The buffer diluent contains low concentrations of toxic sodium azide as a preservative, therefore avoid skin/eye contact and/or ingestion.
- The sample material must be seen as potentially infectious and disposed of accordingly, together with the used test-kit components.

8. TEST PRINCIPLE

The **FASTest® IgG** bovine is designed on latest competitive lateral flow technology based on an immunochromatographic sandwich principle technique.

Due to the competitive test principle, the intensity of **T** is reciprocally proportional to the concentration of bovine IgG in the sample!

Depending on the concentration of bovine IgG antibodies (AB) of the sample (**sample AB**), there is a more or less strong binding of the **detector ABs** (highly specific, polyclonal, mobile, gold labelled) to the **capture ABs** (fixed monoclonal bovine IgG antibodies) on the **TEST line T**.

In case of missing/low AB concentration below the cut-off of the sample, no/few sample AB bind to the detector ABs. These free detector ABs are captured in high amounts by the capture ABs on **T** and produce a more or less weak to strong intensive **T**.

In case of sample AB concentration, many sample AB are bound to the detector ABs. These then cannot bind to **T** anymore. **T** is only weak to not visible.

A correct test procedure will be indicated by a second pink-purple **CONTROL line C**.

Please note the two light-blue quality control lines on **T** and **C**. They indicate that the test membrane is of good quality and can be used.

9. INFORMATION FOR THE INTERPRETATION

- The interpretation of the test result should always be based on anamnestic and clinical data as well as the therapy and prophylaxis possibilities.
- Any non-described colour or contour variation of **T** and **C** (e.g. greyish, shadow-like lines) has to be considered as unspecific reaction. The test result is invalid, the test should be repeated.
- Due to red hemoglobin background, caused by hemolytic blood samples, the visibility of **T**, especially in case of weak positive samples, could be from worse to not visible.

Negative test result:

- Suspicion on partial to complete FPT/lack of IgG (IDS)
- IgG substitution (colostrum administration) required!

Positive test result:

- adequate passiver transfer/no lack of IgG