

Fast + Simple  
Focused on Veterinary Diagnostics



Specialists in Veterinary  
Laboratory Supplies

Exclusive UK distributor  
For MegaCor Diagnostik  
info@vetlabsupplies.co.uk  
Tel: 01798 874567



## SURGICUTT® Vet. H ad us. vet.

### Buccal mucosal bleeding time (BMBT) – standardised hemostasis diagnostics on-site

Standardised screening test for evaluation of primary hemostasis disorders of the dog

#### Emergency diagnostics

- Surface bleedings like petechiae, ecchymoses etc.
- Nose bleeding
- Gastrointestinal hemorrhage
- Intraoperative bleeding bias
- Parasitica infectious diseases like ehrlichiosis etc.


#### Therapeutic monitoring

e. g. aspirin medication

#### Screening test before operations

#### Von Willebrand disease screening test

in races like Dobermann,  
Sheltie, Scotch terrier



- Simple standardised test procedure
- Fast test interpretation after few minutes
- Reliable clinical diagnostics
- Storage at room temperature (15–25 °C)
- Long shelf life
- Compact test box with 10 tests



# SURGICUTT® Vet. H ad us. vet.

Hemostasis is a complex physiological reaction on bleedings of all kind. Coagulation disorders can base on primary and secondary defects or causes.

The buccal mucosal bleeding time (BMBT) is the time that passes by after placement of a defined incision by means of **SURGICUTT® Vet. H**, until hemostasis occurs. The determination of the BMBT is an important coagulation marker for validation of the individual primary state of hemostasis.

With help of BMBT, suspicion of disorders of the primary hemostasis can be confirmed quickly and simply. Characteristics normally are surface bleedings, caused by infectious and/or immunological reasoned thrombocytopenias, thrombocytopathies, vasculopathies and the genetically inherited von Willebrand disease (vWD) type 1–3.

With help of **SURGICUTT® Vet. H**, by a standardised depth and length of the cut, the veterinarian is enabled on-site to obtain constant and reproducible BMBT results and to give a specific statement about the status of hemostasis in the dog and cat.

## Test procedure



**SURGICUTT® Vet. H** for the dog

Incision: **5.0 mm wide and 1.0 mm deep**

1. Remove **SURGICUTT® Vet. H** from blister, do not touch aperture slot
2. Remove safety clamp
3. Fold lip outside and fix safely
4. Place **SURGICUTT® Vet. H** plainly on the mucosa
5. Press the release button of the incision lancet
6. Start the stop watch immediately after beginning of bleeding
7. Allow undisturbed bleeding of cut surface
8. Absorb draining blood every 10 seconds with the filter paper without touching the cut surface
9. Stop the stop watch immediately after bleeding stands and read the buccal mucosal bleeding time (**BMBT**)

## Test interpretation

### BMBT (buccal mucosal bleeding time)



#### 1. Physiological standard

Dog < **4.0 minutes**

#### 2. Prolonged bleeding time

- Von Willebrand disease
- Thrombocytopathy
  - inherited
  - acquired, e.g. by aspirin
- Thrombocytopenia
  - immune mediated
  - infection mediated, e.g. ehrlichiosis
  - bone marrow disease
- Vascular disease

For evaluation of further suspicion of disorders of the secondary hemostasis, the use of **ACT-VETube** is recommended for the determination of the activated coagulation time.

Distribution:

EN 05-2019

