

Fast + Simple  
Focused on Veterinary Diagnostics



Specialists in Veterinary  
Laboratory Supplies

# MYKODERMOASSAY DTM ad us. vet.

# MYKODERMOASSAY TRIO ad us. vet.

Exclusive UK Distributor For  
MegaCor Diagnostik  
[info@vetlabsupplies.co.uk](mailto:info@vetlabsupplies.co.uk)



## DERMATOPHYTOSIS – highly infectious skin disease with zoonotic potential

Special agar for the qualitative detection of veterinary relevant dermatophytes  
in pocket pets, pets and farm animals

**Fast and reliable diagnostics (colour indicator) based on  
cultural dermatophyte detection**

**At clinical suspicion  
(spotted patchy areas of alopecia, often non-pruritic)**

**Early initiation of therapy measures (zoonosis)**

**1 (DTM) or 3 (SAB/DTM/ESA) special culture media**

**Optimised microscopic colony and conidia evaluation**



- Simple test procedure with hair, dandruff, scabs and feathers
- Incubation at 25-32 °C
- Fast test interpretation after only 2-3 days
- Reliable clinical diagnostics
- Storage: **DTM** at 2-8 °C, **TRIO** at 15-25 °C
- Long shelf life
- Compact test box  
**MYKODERMOASSAY DTM** with 12 test vials  
**MYKODERMOASSAY TRIO** with 5 petri dishes and tooth brushes, respectively

# MYKODERMOASSAY DTM ad us. vet.

# MYKODERMOASSAY TRIO ad us. vet.

Dermatophytoses/ringworm belong to the most frequent infectious dermatoses in pocket pets, pets and farm animals, but also in humans (zoonosis).

They are caused by dermatophytes, filamentous fungi using keratin (skin, hair and claws) as carbon source. The clinically most relevant species are *Microsporum* (*M. canis*, *M. gypseum*, *M. persicolor*), *Trichophyton* (*T. verrucosum*, *T. mentagrophytes*) and *Epidermophyton*. Beside age and immunosuppression, familiar, breeding (especially persian cats) and keeping conditions (breeding, animal shelter, hunting dog, multiple species keeping), travelling, lactation (transmission of infection to puppies) as well as e.g. ectoparasite based diseases and debilitated animals play an important role in developing a ringworm disease. Warm and humid climate is an additional trigger.

In case of clinical suspicion of an ongoing dermatophytosis (spotted, patchy areas of alopecia, often non-pruritic), establishment of a mycological culture using dermatophyte specific media is known to be the most reliable technique.

MEGACOR offers two different test formats of special dermatophyte culture media: The **MYKODERMOASSAY DTM**, a classical dermatophyte test media in slant agar vials format and the **MYKODERMOASSAY TRIO**, an innovative petri dish format with a combination of the three most important dermatophyte culture media.

Both test formats ensure a fast evaluation of the clinically suspected diagnosis through colour change of DTM and ESA media, respectively. Furthermore, ESA is an optimised medium for the reliable visual and microscopic differentiation of dermatophytes in practice. SAB/SDA, known as the "classical fungi universal agar", enables a faster growth, compared to DTM and ESA, as well as an optimal double-sided evaluation of the colonies, based on its transparent medium.

## MYKODERMOASSAY TRIO

Culture medium	Colour indicator	Interpretation
<b>SAB/SDA</b> <b>Sabouraud Agar</b> Classical universal fungus medium	<b>No colour change</b> transparent	<b>Growth of all fungus species</b> → perfect optical differentiation from both sides of petri dish
<b>DTM</b> <b>Dermatophyte test medium</b> Specific dermatophyte medium	<b>Colour change from day 2-3 on</b> 	<b>Growth of dermatophytes</b> → optical differentiation of dermatophyte species
<b>ESA</b> <b>Enhanced sporulation agar</b> Specific dermatophyte medium → optimal colony, macro- and microconidia growth	<b>Colour change from day 2-3 on</b> 	<b>Growth of dermatophytes</b> → perfect microscopical differentiation of dermatophyte species

Dermatitis can be caused or intensified by various reasons. Therefore, with skin diseases of unclear aetiology, a diagnostic clarification of a *Malassezia* dermatitis (Pityrosporie) is always recommended.

With the **MYKODERMOASSAY MALASSEZIA** special agar for the qualitative detection of *Malassezia* species in dog or cat, this is achieved easily and reliably in your practice.

Distribution:

UK



Specialists in Veterinary Laboratory Supplies

EN 09-2016

