

The Euromex BioBlue

Binocular Microscope



Specialists in Veterinary
Laboratory Supplies

The Euromex BioBlue Binocular microscope is designed to offer the finest optical components and guaranteed performance, which means clear bright images at every magnification. It is a high-quality build with a 10 year warranty covering mechanical parts and a 2 year warranty on all electrical parts. Of critical importance is its high-quality LED illumination and its semi-plan objective light system. These features allow the system to be used for the vast majority of all veterinary microscopy procedures including high power (1000X) magnifications (e.g. Gram-stained smears & cytology).



MICROSCOPE SPECIFICATIONS

- Head: Binocular
- Nosepiece :Reversed nosepiece for 4 semi plan objectives
- Objectives: Supplied with semi plan 4x/0.10, 10x/0.25, S40x & S100x oil immersion
- WF10x/18mm eyepiece(s)
- Stage: Integrated X-Y mechanical stage
- Focusing: Coaxial coarse and fine adjustments with 200 graduations. Precision 2 μ m per graduation, 0.4 mm per rotation.Total travel approximately 23 mm with rack stop
- Condenser: Height adjustable Abbe condenser with iris diaphragm and filter holder
- Illumination: 1 W NeoLED™ illumination for increased light output
- Power: External mains adapters (100-240V)to 4.5V/1A
- Ergonomic stand
- Warranty 10 years

WARRANTY

10 Year Warranty on Microscope Mechanics 2 Year Warranty on Electrics

The Euromex BioBlue Binocular Microscope SKU 4260

Suitable for most laboratory procedures including those that require higher magnifications. Including worm egg identification, hair/skin scrapes, urine sediments Blood smear stains & manual cell counts.

Visit our website to see our full range of Veterinary Equipment, Diagnostic Test Kits & Consumables...

www.vetlabsupplies.co.uk