

The Euromex BioBlue

LCD Microscope



Specialists in Veterinary
Laboratory Supplies

The Euromex BioBlue with a 7" LCD screen. This microscope has all the features of the regular BioBlue Microscope but is standard equipped with a 7" HD LCD screen. This screen allows live viewing, save images and videos and even has a basic line measurement function included

MICROSCOPE SPECIFICATIONS

- Head: Binocular with Digital LCD Screen
- Nosepiece :Reversed nosepiece for 4 semi plan objectives
- Objectives: Supplied with semi plan 4x/0.10, 10x/0.25, S40x/ Also S60X or S100x oil immersion are also available to purchase separately.
- WF10x/18mm eyepiece(s)
- Stage: Integrated X-Y mechanical stage
- Focusing: Coaxial coarse and fine adjustments with 200 graduations. Precision 2 µm per graduation, 0.4 mm per rotation.Total travel approximately 23 mm with rack stop
- Condenser: Height adjustable Abbe condenser with iris diaphragm and filter holder
- Illumination: 1 W NeoLED™ illumination for increased light output
- Power: External mains adapters (100-240V)to 4.5V/1A
- Ergonomic stand
- Warranty 10 years



LCD SCREEN

- 7" LCD Screen
- Screen resolution x 1080@ 30 fp
- Saves images and videos directly to microSD card
- Image resolution up to 1844 x 1080 in JPEG format
- Video resolution 1280 x 720 AVI format
- Basic measurement function slide available as option
- Power via external mains adapter lead to 5V/1A plug

WARRANTY

10 Year Warranty on Microscope Mechanics 2 Year Warranty on Electrics & LCD Screen

**7 Day Money Back
Guarantee**

The Euromex BioBlue LCD Microscope SKU 4220-LCD

**Suitable for most laboratory
procedures even those that
require higher magnifications.
Including worm egg
identification, hair/skin
scrapes, urine sediments
Blood smear stains & manual
cell counts.**

Visit our website to see our full range
of Veterinary Equipment, Diagnostic
Test Kits & Consumables...

www.vetlabsupplies.co.uk