

The Euromex Microblue

Monocular Microscope



Specialists in Veterinary
Laboratory Supplies

The MicroBlue with a monocular head. This microscope has all the features of the regular MicroBlue Microscope

MICROSCOPE SPECIFICATIONS

- Head: Monocular
- Nosepiece: Quadruple revolving with click stops
- Objectives: Achromat 35mm parfocal 4X, 10X, S40X facing away from operator.
- Stage: Mechanical 115 x 100mm with integrated X/Y travel 55 x 22mm with spring slide holder (accepts McMaster worm egg counting slide).
- Focusing: Coaxial coarse and fine focus controls, 200 graduations, 15 um per graduation, 3mm per rotation, total travel 15mm. Also, friction focus adjustment on course control and adjustable stage stop to protect objectives
- Condenser: Abbe 1.25 condenser, height adjustable with iris diaphragm and filter holder fitted with white filter.
- Illumination: Adjustable 1W LED.
- Power: External mains adapters (100-240V) to 4.5V/1A or Rechargeable Batteries
- Dimensions: 36cm(H) X 15cm(W) X 20cm(D)
- Weight: 2.4 kilos

OPTIONAL EXTRAS (available to purchase separately)

- x20, x60 & x100 Oil Objectives
- 7" LCD Screen
- Binocular Head

WARRANTY

5 Year Warranty on Microscope

**Also Available As Part of Our
Discounted DIY FEC Kits**



The Euromex Microblue Monocular Microscope SKU 1051

**Suitable for laboratory
procedures which do not
require magnifications of
greater than 600X e.g. worm
egg identification, hair/skin
scrapes, urine sediments &
manual cell counts.**

Visit our website to see our full range
of Veterinary Equipment, Diagnostic
Test Kits & Consumables...

www.vetlabsupplies.co.uk