



Hawksley 1400 Centrifuges

Single or Dual Speed Option (75mm microhaematocrit and 2ml microtubes)

The Hawksley Haematospin 1400 Centrifuges are built in the UK and are made from durable and high quality components. These centrifuges are ideal for Veterinary laboratories and practices specifically designed to produce reliable and accurate results on blood and urine. The classic design of these centrifuges is used world-wide with guaranteed support for all servicing and parts.

There is a choice of rotors to allow micro chemistry preparations and it can be specified as either a single or dual speed machine, allowing you to purchase whichever is most suitable for your practice.



Product Information

- ❖ High-speed drive developing RCF 13,000g for fast, accurate, haematocrits, second low speed option (dual speed centrifuge only) developing RCF 3,300g
- ❖ High quality balanced motor for vibration-free, quiet operation without temperature increase
- ❖ Excellent acceleration, auto-braking and smooth running
- ❖ Variable Timer Control, 0-15 minutes, graduated in 1 minute increments
- ❖ Cast aluminium body giving total stability and low noise characteristics
- ❖ Simple intuitive operation suitable for any level of operator
- ❖ Range of easy to change and versatile rotors
- ❖ Characteristic classic styling
- ❖ Lockable steel lid and rugged components that will withstand years of hard service
- ❖ **Dual Speed Option:** Fast – 13,000g for blood Slow – 2,200g for urine

Technical Specifications

Speed	11,800rpm with 24-place rotor
RCF	14,000g with 24-place rotor
Acceleration	0-11,800rpm in 60 seconds, with 24-place rotor
Brake	Automatic
Noise	62dB running with 24 place rotor
Motor	High precision AC/DC series wound motor, brush gear, computer wound coil winding, balance armature and high speed bearings
Timer	0-15 minute timer, graduated in 1 minute increments. Automatic switch at zero
Weight	5.1kg net (without rotor)
Dimensions	23.5cm height x 27cm diameter, small footprint 16cm diameter
Construction	Aluminium casing, steel lid
Finish	Body cream white baked enamel, lid matt black baked enamel







Price List

Cat No	Product Description	Price
1513	Hawksley Vetspin 1400 Dual Speed Centrifuge with Multi-Combi Rotor	£965.00

Cat No	Product Description	Price
1500	Hawksley Haematospin 1400 Single Speed Centrifuge without Rotor	£690.00
1507	Hawksley Haematospin 1400 Single Speed Centrifuge + 16 place Combi Rotor	£815.00
1505	Hawksley Haematospin 1400 Single Speed Centrifuge + 8 place Acetal Rotor	£768.00
1506	Hawksley Haematospin 1400 Single Speed Centrifuge + 24 place Mhct Rotor	£815.00

Rotors

All the above rotors are available separately, details and prices below.

	24 Place Haematocrit Rotor for 75mm x 59µl capillary tubes, in precision-machined aluminium Cat No 1504 Price £128.82
	Multi-Combi Rotor for 8 Haematocrit, plus 16 micro tubes in 3 sizes, in precision-machined aluminium Cat No 1600 Price £128.82
	Combi Rotor for 16 Haematocrit, plus 16 1.5/2ml micro tubes, in precision aluminium Cat No 1503 Price £128.82
	8 place Acetal Rotor for 1.5/2ml micro tubes Cat No 1554 Price £77.66

Prices subject to VAT

**We also have a specialist Hawksley Service & Repair department
Plus...**

**We also sell centrifuge laboratory consumables & much more!
Contact us for a consumables price list.**

Aug 2015