



The Vetlab Premiere 125 Binocular Laboratory Microscope

The Vetlab Premiere 125 Binocular Microscope has been specially selected for its solid robust build and its reliable performance. This is a classically designed, basic compound microscope that will stand the test of time in any laboratory.

Features:

- **Optical System:** Finite 160mm.
- **Viewing Head:** Sliding Binocular (Jentsch). Interpupillary distance 55mm – 75mm inclined at 45° and rotatable by 360°. Diopters on each eyepiece aperture.
- **Eyepieces:** One Set of Wide Field (WF) 10x/ 18mm eyepieces.
- **Nosepiece:** Quadruple.
- **Objectives:** Finite Achromat 4X, 10X, 40X (S), 100X (S) (O). Parfocal and mounted towards the operator (convenient for cleaning lenses).
- **Focus Controls:** Coaxial coarse and fine focus controls, range 30mm; fine focus interval 0.002mm. Focus/rack stop control and Focus Tension control.
- **Stage:** Double layer mechanical stage 140mm x 140mm, moving area 75 x 45mm with a spring clip single slide holder.
- **Condenser:** Abbe Condenser N.A. 1.25. Iris diaphragm and filter holder. Rack and pinion height adjustment.
- **Light Source:** Halogen lamp 6V/20W. On/Off and variable intensity switches.
- **Starter Pack:** User Manual, mains power lead, spare halogen bulb, spare fuse, microscope slides, coverslips, lens tissues, immersion oil and dust cover.
- **Optional Extras:** 60X Dry Achromat Objective.



3-Year Warranty CE Certified

AVAILABLE ON 7 DAY SALE OR RETURN – UK MAINLAND ONLY- TERMS & CONDITIONS APPLY

Cat No	Description	Price
1932	Vetlab Premiere 125 Binocular Microscope with four Achromat Objectives	£285.00

Price subject to VAT – FREE DELIVERY UK MAINLAND