



The Vetlab Premiere 2000 Binocular & Trinocular LED Laboratory Microscope

The Vetlab Premiere 2000 LED Veterinary Microscope is our top of the range microscope designed with the professional microscopist in mind. Its superb modern, ergonomic design makes it very comfortable and easy to use. The attention to detail used in the design and construction of this microscope, will ensure a high level of comfort during long periods of microscopy examinations.

Features:

- **Optical System:** INFINITY
- **Viewing Head:** Binocular or Trinocular (Seidentopf) fitted with Infinity tube lens. Interpupillary distance 48mm – 75mm, inclined 30° rotatable 360°. Diopters on each eyepiece aperture.
- **Eyepieces:** One set of **Extra Wide Field (WF) 10X/22mm** eyepieces.
- **Nosepiece:** Quadruple.
- **Objectives:** High contrast, chromatic-free full **PLAN INFINITY** objectives 4X, 10X, 40X (S), 100X (S)(O). Parfocal and mounted away from the operator (convenient for changing slides).
- **Focus Controls:** Coarse and fine controls, range 38mm; fine focus division 0.001mm. Focus/rack stop control and Focus Tension control.
- **Stage:** Double layer large mechanical stage 216mm x 150mm, moving area 75mm x 55mm with a spring clip single slide holder. Low positioned X-Y controls.
- **Condenser:** Abbe condenser N.A. 1.25, Iris Diaphragm and filter holder, Centring adjustment screws. Rack and pinion height adjustment.
- **Light Source:** Built-in LED 3W with Field Diaphragm for Kohler Illumination.
- **Starter Pack:** User Manual, mains power lead, spare fuse, microscope slides, coverslips, lens tissues, immersion oil and dust cover.
- **Optional Extras:** 60X Dry Plan Infinity Objective.
Trinocular Head with eyepiece and C-Mount.
Dark ground illumination.
Phase contrast.



3-Year Warranty **CE Certified**

AVAILABLE ON 7 DAY SALE OR RETURN – UK MAINLAND ONLY – TERMS & CONDITIONS APPLY

Cat No	Description	Price
1946	Vetlab Premiere 2000 LED Binocular Microscope	£599.00
1947	Vetlab Premiere 2000 LED Trinocular Microscope	£640.00

Price subject to VAT – FREE DELIVERY UK MAINLAND