

# The Vetlab Premiere 250

## LED Binocular Microscope



Specialists in Veterinary  
Laboratory Supplies

The Vetlab Premiere 250 veterinary LED microscope, has a superb ergonomic design and innovative features. These features make it the perfect choice for routine diagnostic examinations in all veterinary and medical laboratories as well as pharmaceutical research. The Premiere 250 professional microscope is supplied with a binocular head.

### Features:

#### Optical System:

Finite 160mm.

#### Viewing Head:

Binocular (Siedentopf). Interpupillary distance 48mm – 75 mm, inclined 30°, rotatable 360°. Diopters on each eyepiece aperture.

#### Eyepieces:

One Set of Wide Field (WF) 10X/18mm eyepieces.

#### Nosepiece:

Quadruple.

#### Objectives:

Finite Achromat 4X, 10X, 40X (S), 100X (S)(O). Parfocal and mounted towards operator (convenient for cleaning lenses).

#### Focus Controls:

Coaxial coarse and fine focus controls, range 36mm; fine focus division 0.002mm. Focus/rack stop control and Focus Tension control.

#### Stage:

Double layer mechanical stage 142mm x 132mm Moving area 75mm x 50mm with a spring clip single slide holder.

#### Condenser:

Abbe condenser N.A. 1.25. Iris diaphragm and filter holder. Centring adjustment screws. Rack and pinion height adjustment.

#### Light Source:

Built-in, LED 3W. On/Off and variable intensity switches.

#### Included:

User Manual, mains power lead, spare fuse.

#### Free:

Microscope slides, coverslips, lens tissues, immersion oil and dust cover.

#### Optional Extras:

60X Dry Achromat Objective.



### **The Vetlab Premiere 250 LED Binocular Microscope**

3 Year Warranty

Free Delivery : UK Mainland

Free Microscope Consumables

Available On 7 Day Sale Or Return

(See Terms & Conditions)

CE Certified

Visit our website to see our full range of Veterinary Equipment, Diagnostic Test Kits & Consumables...

- Special Offers
- New Products

[www.vetlabsupplies.co.uk](http://www.vetlabsupplies.co.uk)