

The Vetlab Premiere 260

LED Binocular Microscope



Specialists in Veterinary
Laboratory Supplies

The Vetlab Premiere 260 LED microscope offers exceptional value for money with its robust build, modern ergonomic features and high specifications. It is ideal for the small to medium size practice laboratory where a variety of procedures are required on a frequent basis. It has a powerful built-in 3W LED light source and can also run on batteries for cordless/field use. A Trinocular version is available for photomicrography/teaching applications.

Features:

Optical System: FINITE

Viewing Head:

Binocular (Siedentopf)

Interpupillary distance 55mm – 75mm, inclined 30° and rotatable 360°

Diopters on each eyepiece aperture.

Eyepieces:

One set of Wide Field (WF) 10X/18mm eyepieces.

Nosepiece:

Quadruple.

Objectives:

ACHROMAT objectives 4X, 10X, 40X (S), 100X (S) (O).

Parfocal and mounted away from operator (convenient for changing slides).

Focus Controls:

Ergonomically designed coaxial coarse and fine controls, range 25mm; fine focus division 0.002mm. Focus/rack stop control and Focus tension control.

Stage:

Double layer mechanical stage 140 x 130mm, moving area 75mm x 50mm with a spring clip single slide holder. Ergonomically designed, low positioned X -Y controls.

Condenser:

Bright Field Abbe Condenser N.A. 1.25, Iris Diaphragm and filter holder. Centring adjustment screws. Rack and pinion height adjustment.

Light Source:

Built-in LED 3W 220V, power indicator light. On/Off and variable intensity switches. **Battery compartment for 3 x AA batteries for cordless/field use.**

Included:

User Manual, mains power lead, spare fuse.

Free:

Microscope slides, coverslips, lens tissues, immersion oil and dust cover.

Optional Extras: 60X Dry Achromat Objective, Plan 10X/20mm WF Eyepieces with rubber eye-guards, Heated Stage (for semen viewing), Photographic Accessories, Soft Sided Carrying Case.



The Vetlab Premiere 260

LED Binocular Microscope

Cat. No. 1964

2 Year Warranty

Free Microscope Consumables

CE Certified

Vetlab Supplies Ltd Unit 13 Broomers Hill Park Broomers Hill Lane Pulborough West Sussex RH20 2RY

email info@vetlabsupplies.co.uk www.vetlabsupplies.co.uk Telephone 01798 874567 Fax 01798 874787