

# The Vetlab SP150 Plan Infinity Trinocular Microscope



Specialists in Veterinary  
Laboratory Supplies

The Vetlab SP150 is a Plan Infinity Trinocular Microscope with a high level specification at a very reasonable price. It has a modern, ergonomic design and the Plan Infinity optical system combined with the bright LED Kohler principal illumination makes it ideal for photomicrography and phase contrast applications. Accessories include a pilot style carrying case.

## Features:

### Optical System:

INFINITY

### Viewing Head:

Trinocular (Seidentopf) fitted with Infinity tube lens. Inclined 30°, rotatable 360°.

### Eyepiece

Wide Field (WF) 10X/18mm eyepieces with dioptic and interpupillary adjustment.

### Nosepiece:

Quintuple.

### Objectives:

Infinity Planachromatic DIN standard objectives 4X, 10X, 20X, 40X (S), 100X (S) (O).

### Focus Controls:

Low position coaxial coarse and fine focus controls with tension adjustment on the focus mechanism. User defined Focus Safety Stop.

### Stage:

Large mechanical stage with coaxial drop down controls and double Vernier scales, spring clip single slide holder.

### Condenser:

Full sub-stage Abbe Condenser (N.A. 1.2), Iris Diaphragm with full lateral adjustment, filter carrier and full rack and pinion height adjustment.

### Light Source:

Built-in LED light source with Field Diaphragm for Kohler Illumination. Rheostat controlled.

### Included:

User Manual, mains power lead, spare fuse.

### Optional Extras:

60x Plan Objective, CCTV, Measuring and Photography Adapters  
Phase Contrast Set, Dark Field Condenser, Pilot Style Carrying Case  
Other Accessories available



## The Vetlab SP150 Trinocular Microscope

2 Year Warranty

CE Certified

Visit our website to see our full range of Veterinary Equipment, Diagnostic Test Kits & Consumables...

- Special Offers
- New Products

[www.vetlabsupplies.co.uk](http://www.vetlabsupplies.co.uk)