

The Vetlab SP22/22D

LED Monocular Microscope



Specialists in Veterinary
Laboratory Supplies

The Vetlab SP22 LED Monocular Microscope is an ideal starter microscope suitable for the hobby user or amateur naturalist. It has many features normally found on more expensive models, representing fantastic value for money. This microscope is also available with an inbuilt digital camera (SP22D).

Please note neither of these microscopes will accommodate a McMaster Slide.

Features:

Optical System:

FINITE

Viewing Head:

Monocular, rotatable 360°

Eyepiece:

Wide Field (WF) 10X/18mm eyepieces.

Nosepiece: Quadruple.

Objectives:

Finite Achromatic RMS standard objectives 4X, 10X, 40X (S), 100X (S) (O).

Focus Controls:

Coaxial coarse and fine focus controls. Tension adjustment on the focus mechanism. Focus Safety Stop.

Stage:

Mechanical stage with a spring clip single slide holder. Coaxial drop down controls and double Vernier scales.

Condenser:

Full sub-stage Abbe Condenser (N.A. 1.2) with Iris Diaphragm and rack and pinion height adjustment.

Light Source:


Built-in LED light source. Rheostat controlled.

Included: User Manual, mains power lead, spare fuse.

Optional Extras: Aluminium Carrying Case

SP22D

This microscope has the same optical performance and mechanical build as the SP22 but incorporated in the head is a 1.3M pixel (1280 x 1024) CMOS chip digital camera. It shares the image with the eyepiece and is attached to a laptop or PC via a USB lead (included). The supplied Scopeimage software is suitable for Windows 7, 8 and 10 but not for Apple OS. Step by step instructions provided.



The Vetlab SP22/22D Monocular Microscope

Cat. No. 1953 SP22

Cat. No. 1955 SP22D

**Cat. No. 1956 Carrying Case
(SP22 only)**

2 Year Warranty

CE Certified

Visit our website to see our full range of Veterinary Equipment, Diagnostic Test Kits & Consumables...

- Special Offers
- New Products

www.vetlabsupplies.co.uk