

The Wedmore SP03

Monocular Microscope



Specialists in Veterinary
Laboratory Supplies

The Wedmore SP03 is a low cost Plano/Concave Mirror Monocular Microscope with a robust design ideal for use in education establishments and for the dedicated amateur naturalist. It has traditional mirror illumination for those with no access to electricity or who need a portable option and comes with its own carrying case, making it even better value for money.

Features:

Optical System:

FINITE

Viewing Head:

Monocular

Eyepiece

Wide Field (WF) 10X/18mm eyepieces.

Nosepiece:

Quadruple.

Objectives:

Achromatic DIN standard objectives 4X, 10X, 40X (S), 100X (S) (O).

Focus Controls:

Coaxial coarse and fine focus controls. Tension adjustment on the focus mechanism. Focus Safety Stop.

Stage:

Mechanical stage with a spring clip single slide holder. Coaxial drop down controls and double Vernier scales.

Condenser:

Full sub-stage Abbe Condenser (N.A. 1.2) with Iris Diaphragm and filter holder. Rack and pinion height adjustment.

Light Source:

Plano/Concave Mirror Illumination.

Included: User Manual, mains power lead, spare fuse.



The Wedmore SP03 Monocular Microscope

Cat No: 1951

2 Year Warranty

CE Certified

Visit our website to see our full range of Veterinary Equipment, Diagnostic Test Kits & Consumables...

- Special Offers
- New Products

www.vetlabsupplies.co.uk

Comes complete with Carrying Case

Vetlab Supplies Ltd Unit 13 Broomers Hill Park Broomers Hill Lane Pulborough West Sussex RH20 2RY
email info@vetlabsupplies.co.uk www.vetlabsupplies.co.uk Telephone 01798 874567 Fax 01798 874787