

The Wedmore SP14

Monocular Microscope



Specialists in Veterinary
Laboratory Supplies

The Wedmore SP14 is a low cost Monocular Microscope with an LED brightness controlled illumination system. It has a robust design ideal for use in education establishments and for the dedicated amateur naturalist. An excellent microscope without equal in its price range, it comes with its own carrying case which makes it even better value for money.

Features:

Optical System:

FINITE

Viewing Head:

Monocular

Eyepiece

Wide Field (WF) 10X/18mm eyepieces.

Nosepiece:

Quadruple.

Objectives:

Achromatic DIN standard objectives 4X, 10X, 40X (S), 100X (S) (O).

Focus Controls:

Coaxial coarse and fine focus controls. Tension adjustment on the focus mechanism. Focus Safety Stop.

Stage:

Mechanical stage with a spring clip single slide holder. Coaxial drop down controls and double Vernier scales.

Condenser:

Full sub-stage Abbe Condenser (N.A. 1.2) with Iris Diaphragm and filter holder. Rack and pinion height adjustment.

Light Source:

Built-in LED light source. Rheostat controlled.

Included: User Manual, mains power lead, spare fuse.



The Wedmore SP14 Monocular Microscope

Cat No: 1952

2 Year Warranty

CE Certified

Visit our website to see our full range of Veterinary Equipment, Diagnostic Test Kits & Consumables...

- Special Offers
- New Products

www.vetlabsupplies.co.uk

Comes complete with Carrying Case

Vetlab Supplies Ltd Unit 13 Broomers Hill Park Broomers Hill Lane Pulborough West Sussex RH20 2RY
email info@vetlabsupplies.co.uk www.vetlabsupplies.co.uk Telephone 01798 874567 Fax 01798 874787